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WINCARE

**User manual**  
**Relief devices and foam positioning system**  
**POSIMO**

**DEVICE INTENDED USE**

Disability prevention and compensation.

**GENERAL INDICATIONS**

Bedridden patients with loss of autonomy and mobility with impairment of functional or structural integrity. Helps in the prevention and treatment of pressure ulcers.

**TARGET POPULATION**

Adults and adolescents 12 years and older

**CONTRAINDICATIONS**

Severe deterioration of cardio-respiratory status (VCDA/ 50). Acute and/or procedural pain. Unstabilised fractures on medical advice.

**ADVERSE SIDE EFFECTS**

Excessive perspiration on direct skin contact with the cover. Position-related discomfort. Acute and procedural pain.

**Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the device poses a serious risk or has been tampered with.**

**SPECIFIC (CONTRA)INDICATIONS / MEDICAL DEVICE COMPOSITION / TECHNICAL DATA**

Reference	Indications	Contraindications	Foam	Cover
VCDA/50	Preserve the functional (respiratory, analgesic) and structural (especially cutaneous) integrity of the pelvis and thorax during prolonged bed rest.	Adverse clinical condition to mobilization (post-traumatic), pain related to support. Agitation.	PU VE	58% Polyester / 35% Polyurethane- (If CIC Polycarbonate with Ag <sup>+</sup> ions) Oeko-Tex 100 class I
VTLA/HIPH, VTLAXL/HIPH	To preserve the functional integrity (notably analgesia) and structural integrity (notably cutaneous) of the heel in the patient without risk of equinus. Lower limb arteritis.	-	PU VE	
VDT41A/CIC, VDT45A/CIC	To preserve the functional integrity (notably analgesia) and structural integrity (notably cutaneous) of the heel in the patient without risk of equinus.	-	PU VE	
VTAC410/CIC, VTAC450/CIC	To preserve the functional integrity (notably analgesia) and structural integrity (notably cutaneous) of the heel in the patient with risk of equinus.	Malleolus arterial ulcer. Severe lower limb oedema.	PU HR	
VTM54	Abduction position of the lower limbs when turning over during changing or certain nursing or medical procedures.	Agitation. Severe lower limb oedema.	PU HR	

\*Consult the terms and conditions for handling products in the EU member state whether they are used alone or in combination.

**CLINICAL BENEFIT, DEVICE PERFORMANCE, MECHANISM OF ACTION**

CHARACTERISTICS IN TERMS OF DEVICE PERFORMANCE: high resilience (HR) or viscoelastic (VE) polyurethane (PU) foam contained in a protection that withstands penetration of water, water vapour, that is bi-elastic that immerses those parts in contact and relieves or eliminates body pressure points that risk causing bedsores in the lying position.

EXPECTED CLINICAL BENEFITS: Maintaining a sustainable position. Maintenance of tissue oxygenation by relieving or removing pressure from body areas vulnerable to contact with the mattress.

INFORMATION FOR HEALTHCARE PROFESSIONALS: Monitor the condition of the patient's skin in contact with the support several times a day. Use pressure relief devices or positioning systems in patients with pressure ulcer(s). Change position at least every 8 hours on a motorized alternating pressure air support, every 4 hours on a viscoelastic foam mattress and every 2 to 3 hours on all others.

**PREREQUISITES BEFORE USE AND INSTRUCTIONS FOR USE**

TRAINING AND QUALIFICATION OF THE DEVICE USER: user training in the positioning and installation of patients must be carried out by persons trained and approved by the economic operators concerned (including safety and non-compliance reporting).

DEVICE installation: The technical assistance is delivered with its removable zipped cover (except VTMA54). The products are ready to install. The surface of the high resilience foam (VTMA54) or viscoelastic foam (others) must be positioned in contact with the patient's body. The cover must be dry. Preventive maintenance: Carry out a regular visual check of the condition of the foam: visible sagging of material and a non-homogeneous and very slow return of the foam are ageing criteria that compromise the properties of the support. Annual inspection of the condition of the cover (surface aspect and zippers) by exposing the inner surface to a light source: no holes and/or tears. Replace the cover if its surface appearance changes.

**CLEANING AND DISINFECTION**

Cover						
Moderate washing up to 90°C	Maximum chlorine concentration of 5000ppm	Do not iron	Do not dry clean	Dry with reduced thermal restriction	Use an approved surface disinfectant/detergent	

**INFORMATION ON APPROPRIATE PROCESSES TO ALLOW ITS REUSE**

The product must exhibit physical and bacteriological cleanliness. Item treated with a biocide substance at no risk to users.

Avoid wet contact with the foam. Prohibit all scouring, stripping products or solvents and "piercing-sharp" objects in direct contact with the cover.

**WARNINGS, PRECAUTIONS FOR USE, MEASURES REQUIRED****PRECAUTIONS FOR USE**

Manually check spine/pelvis alignment and trochanter (VCDA / 50) pressure removal. The bed surface is horizontal.

In the case of an established pressure sore (from stage I), avoid direct contact with the device. If the patient is agitated / confused and / or if the heel relief is not followed, use VTLA/H or VTAXL/HIPH.

**WARNINGS**

At each change in position, check the appearance of the skin in direct contact with the system surface (insides of the knees, heels, malleoli if external rotation of a lower limb, etc.). At each 30° semi-lateral position change, check for redness at the foot outer edge and the malleolus in contact with the mattress, of the ear in contact with the technical assistance. Excessive perspiration of the zones in extended direct contact with the cover surface (towel interface).

**ACTION REQUIRED**

- Read the recommendations of the basic or therapeutic mattress manufacturer on the use of these relief and associated positioning systems.
  - When applicable, check that the bed railing height is compatible with the mattress thickness or mattress and over-mattress thickness (at least 22 cm between the top of the mattress and the top of the railings) and that the positioned body parts are not in contact with the medical bed structure and its accessories;
  - Important: changing patient position is not sufficient to prevent bedsores; other preventive measures are also essential:
    - use of a therapeutic mattress adapted to the patient risk, mobility, activity;
    - maintaining skin hygiene and avoiding maceration;
    - in the event of incontinence, changing protection regularly;
    - examine skin condition (or have it examined) every day;
    - make sure food intake is sufficient and drink regularly and in sufficient quantities.
- If any of these measures cannot be followed, your doctor or nurse must be informed as quickly as possible.
- It is important to limit over-thickness between the body and the support as much as possible, except for a comfort cover compatible with positioning systems, body clothing and a potential complete change. Favour loose cotton clothing and if possible without stitching in the support zone.
  - Make sure there are no foreign objects such as: tubes, crumbs, fat, etc. before positioning the person,
  - Important: consult your prescriber within one month of purchasing these relief or positioning systems.

**CIRCUMSTANCES IN WHICH USERS SHOULD CONSULT WITH A HEALTHCARE PROFESSIONAL**

Immediately notify your doctor or nurse of all abnormal events such as fever, pain or redness or whitening of pressure points (head, shoulder, back, hip, shoulder blade, pelvis, heel, etc.).

**INFORMATION ON ALL KNOWN RESTRICTIONS ON THE COMBINATION WITH SYSTEMS AND EQUIPMENT**

The use of protection other than that provided and specified by WINNCARE may result in the system malfunctioning and/or render it incompatible with safe use and the preservation of its performances.

**STORAGE, HANDLING, DISPOSAL****STORAGE AND USAGE CONDITIONS**

POSIMO products should preferably be stored flat, away from direct light and excessive humidity.

	Use	Storage
<b>Temperature range</b>	+15°C to +45°C +59°F to +113°F	-5°C – 60°C +41°F to +140°F
<b>Humidity range</b>	30% - 70%	10% - 90%

**SERVICE LIFE**

Relief devices or positioning systems have an estimated service life of 3 years.

**PRODUCT DISPOSAL**

Do not dispose of the product outside of the dedicated area. Respect the established recycling systems in your country.

**GUARANTEE**

The guarantee duration for the POSIMO range is 2 years. This guarantee starts from the date of purchase of the product from your distributor. This guarantee does not cover normal wear and tear of the product and its protective cover and is not a substitute for legal guarantees. Please contact your distributor about a defective product, specifying the batch number of the product. They will carry out the necessary procedures with WinnCare to proceed with either repair or exchange.

Quality management system: ISO 13485:2016 - Certifications covering QHE (ISO 9001: 2015 – ISO 14001:2015 – OHSAS 18001: 2017)

**Class I medical device according to Regulation (EU) 2017/745.****Contents of package**

1 relief or positioning system  
 1 user manual

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Siotec AG, Sägestrasse 75, CH-3098 Köniz