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**DEVICE INTENDED USE**

Pressure sore prevention in prolonged sitting position.

**INDICATIONS**

For bedridden patients with mobility impairment with "Low to Medium" risk of pressure ulcers (according to the Braden scale or any other approved scale and in accordance with medical opinion).

**CONTRAINDICATIONS**

Patient with sagittal or frontal stability disorders. Paralysis of the trunk and/or lower limbs. Paralysis of the trunk and/or lower limbs, history of ischial or sacrococcygeal pressure ulcer, lower limb tenderness, spasticity, vascular amputation. Patient weight level below or above the limits defined in the technical data table.

**ADVERSE SIDE EFFECTS**

Trunk dynamic instability. Pelvic slippage in the sagittal plane.

**Any serious incident occurring in connection with the device must be notified to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. Inform the competent authority if you consider or have reason to believe that the device poses a serious risk or has been tampered with.**

**MEDICAL DEVICE COMPOSITION/TECHNICAL DATA**

Models	VAB (GELTOP)				VGLS (GELSCAR)	
	VAB4040/PUHD	VAB4345/PUHD	VAB4040/CIC	VAB4345/CIC	VGLS01C	VGLS01C/4345
Sizes	40x40x5,5 cm	43x45x5,5 cm	40x40x5,5 cm	43x45x5,5 cm	40x40x3,5 cm	43x45x3,5 cm
Specific gel	Viscofluid Coatex DV1163 Thickness 1 cm	Viscofluid Coatex DV1163 Thickness 1 cm	Viscofluid Coatex DV1163 Thickness 1 cm	Viscofluid Coatex DV1163 Thickness 1 cm	Viscofluid Coatex DV1163	Viscofluid Coatex DV1163
Specific Lower Foam	PU HR 34 kg/m <sup>3</sup>	PU HR 34 kg/m <sup>3</sup>	PU HR 34 kg/m <sup>3</sup>	PU HR 34 kg/m <sup>3</sup>	-	-
Protection	Removable cover PU51%/Polyester49% Non-slip lower surface	Removable cover PU51%/Polyester49% Non-slip lower surface	Removable cover PU42%/Polyester58% Non-slip lower surface	Removable cover PU42%/Polyester58% Non-slip lower surface	Fixed neoprene cover	Fixed neoprene cover
Patient weight min and max	40 to 95 kg	40 to 95 kg	40 to 95 kg	40 to 95 kg	50 to 110 kg	50 to 110 kg

The references mentioned can be completed with a suffix to constitute a commercial reference (example VA40TH)

**CLINICAL BENEFIT, DEVICE PERFORMANCE, MECHANISM OF ACTION**

**DEVICE PERFORMANCE CHARACTERISTICS:** The viscofluid gel principle as an interface with the seated patient reduces shear force in particular in a dynamic position. The combination of viscofluid gel with a foam base allows relative immersion of the pelvis and the areas opposite bony protuberances (ischium, coccyx, trochanter) to be enveloped.

**PROJECTED CLINICAL BENEFITS:** Maintenance of tissue oxygenation in the cutaneous and subcutaneous tissues in contact with the support.

**INFORMATION FOR HEALTHCARE PROFESSIONALS:** Monitor the condition of the patient's skin in contact with the cushion several times a day. Do not seat patients with ischial and coccygeal pressure ulcers. Check that the seat of the chair is in good condition. In the case of a wheelchair, check the height adjustment of the armrests, footrests, (...) and performance in terms of propulsion.




**PREREQUISITES BEFORE USE AND INSTRUCTIONS FOR USE**

**TRAINING AND QUALIFICATION OF THE DEVICE USER:** user training must be carried out by individuals trained and approved by the economic operators concerned, notably in terms of safety and non-compliance reporting.

**DEVICE installation:** The cushion is delivered with its cover fitted and dry. The product is ready to be installed. The upper surface of the cushion in contact with the patient is made up of viscofluid gel. Attach the Gelscar cushion straps to the chair uprights.

**Preventive maintenance:** Carry out a regular visual check of the cover (surface appearance - absence of holes and/or tears - and zippers) or envelope condition. The presence of a gel leak renders the product unfit for use. For the Geltop: Replace the protective cover if its surface appearance changes.

**CLEANING AND DISINFECTION**

Procedure						
Cover GELTOP	Moderate washing up to 90°C	Maximum chlorine concentration of 5000ppm	Do not iron	Do not dry clean	Dry with reduced thermal restriction	Use an approved surface disinfectant/detergent
Gelscar Envelope	Surface detergent				No	

**INFORMATION ON APPROPRIATE PROCESSES TO ALLOW ITS REUSE**

Before re-using for another patient, the product must exhibit physical and bacteriological cleanliness after application of the above procedures. Do not use the Geltop cushion without a cover. Item treated with a biocide substance at no risk to users. [www.winnicare.fr](http://www.winnicare.fr)

**WARNINGS, PRECAUTIONS FOR USE, MEASURES REQUIRED****PRECAUTIONS FOR USE**

Non-stabilized bone and/or muscle injuries in contact with the support. COPD-type bronchopulmonary disorders with severe reduction in respiratory muscle tone.

A bedsore is a lesion of varying depth in the skin, associated with excessive and prolonged compression of tissue between the body and a support surface. This excessive pressure can prevent blood circulation and lead to pressure ulcers, a dangerous lesion because it develops from the inside of the body to the outside, and is therefore initially invisible. A pressure sore can have several forms: a simple persistent rash for more than a day, hardening of the skin, a wound of varying depth that can in severe cases affect muscles or underlying bone.

Pressure (max) in mmHg before ageing, at the maximum validated weight, measured flat (0°)

SUITABILITY FOR USE	Max weight User	Head	Ischions	Heel
Geltop cushion (mm of Hg)	95 kg	NA	88.5	NA
Gelscar cushion (mm of Hg)	110 kg	NA	82	NA

**STORAGE, HANDLING, DISPOSAL****STORAGE AND USAGE CONDITIONS**

The supports must ideally be stored flat, away from direct light and excessive humidity.

	Use	Storage
Temperature range	+15°C to +45°C / +59°F to +113°F	-25°C to 70°C / +13°F to +158°F
Humidity range	30% - 70%	30% - 95%
Atmospheric pressure	50kPa - 106 kPa // 7.25 Psi / 15.37 Psi	

**SERVICE LIFE**

The cushions have an estimated service life of 6 years.  
(See date on package label).

**PRODUCT DISPOSAL**

Do not dispose of the product outside of the dedicated area. Respect the established recycling systems in your country.

**GUARANTEE**

The guarantee duration for cushions is 2 years. This guarantee starts from the date of purchase of the product from your distributor. In the case of single patient use, the guarantee relates to the support and its protective cover. This guarantee does not cover normal wear and tear of the product and its protective cover and is not a substitute for legal guarantees. Please contact your distributor about a defective product. They will carry out the necessary procedures with our company to proceed with either repair or exchange.

**Technical compliance certificates**

FCBA - N° ESC 17-004	FCBA - N° ESC 19-001
Geltop Cushion	Gelscar cushion

**Class 1 medical device according to Regulation (EU) 2017/745.**

**Management system: Quality: ISO 13485:2016 - Environmental: ISO 14001:2015**

**Contents of package**

- 1 cushion.
- 2 covers for the Geltop cushion. No cover required for the Gelscar cousin.
- 1 user manual

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Processo						
Cover Gelpop	Moderate wash at 90°C	Maximum authorised chlorine concentration 5000ppm	Non stirare	Do not dry clean	Drying with reduced thermal constraints	Use of an authorised surface detergent/disinfectant product
Involucro Gelscar	Detergente per superfici				No	

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(max) Pressure in mmHg before ageing, using the maximum validated weight, measured flat (0°)

IDONEITÀ ALLO SCOPO	Max. weight User	Head	Ischions	Heels
Cuscino in gelpop (mm Hg)	95 kg	NA	88,5	NA
Cuscino in gelscar (mm Hg)	110 kg	NA	82	NA

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This warranty does not cover normal wear and tear and is not a substitution of legal guarantees. Please contact your distributor and present the incriminating product. The distributor will carry out the necessary procedure with our company in order to either repair or replace the product.

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Siotec AG, Sägestrasse 75, CH-3098 Köniz