
# HEBERGEMENT

30 < 130 kg  
67 < 287 lb

70 - 140 cm

**WINNCARE FRANCE, Miniac**  
ZONE ACTIPOLE  
2 RUE DE SAINT-COULBAN  
35540 MINIAC-MORVAN - FRANCE

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### Indication

Device indicated for patients who are bed-bound for a few days, who are able to move alone and without any problems, and who are at not at risk of developing a pressure sore (as per clinical judgement and scales).

### Contraindications

- > Use on medical imaging table and stretcher,
- > Patients whose weight is higher than that which the mattress can support.

### Patient and user target groups

Adult with temporary or permanent loss of independence, or whose height is greater than or equal to 146 cm.

### Indicate the undesired side effects

All serious incidents which are related to the device must be notified to the manufacturer and the competent authority of the member state in which the user and/or patient resides.

Inform the competent authority if you believe, or if you have reason to believe, that the device presents a serious risk, or that it has been tampered with.

## Expected clinical benefits

Contributes to the activities required for life support.

## Information for health professionals

Observe the condition of the patient's skin which is in contact with the support structure several times a day. Check that the bed-base and mattress foam cover are in good condition.

## 4 INSTRUCTIONS FOR USE

### Training and qualifications required for using the device

User training must be carried out by persons who have been trained and approved by the relevant economic operators, especially in terms of safety and reporting of non-conformities.


### Device installation

The mattress comes with a removable zipped cover. The product is ready to be installed. The cover must be dry.

For mattresses which have been delivered in a rolled-compressed condition, you must wait for the decompression time, which can vary from 30 minutes to 2 hours, to pass.


## STEP

Foot/head

Place the mattress on the bed. Lie the patient down.

## 2 COMPOSITION OF THE MEDICAL DEVICE

- ❶ High resilience polyurethane foam mattress (density  $\geq 27$  kg/m<sup>3</sup>).
- ❷ 3-sided zipped cover.
- ❸ Identification label on the outside of the cover.
- ❹ Instructions for use (paper or electronic).

## 3 CLINICAL BENEFIT, PERFORMANCE, MECHANISM OF ACTION

### Performance characteristics of the device

This mattress adapts itself to the different bed-base positions required for the provision of care, the performance of medical procedures, and/or for the comfort of bed-bound patients and basic daily activities (body hygiene, transfers, etc.).

### Preventive maintenance

- > Annually carry out a visual inspection of the state of the foam: the presence of visible subsidence of the material, and an uneven and very slow return of the foam are signs of ageing which compromise the properties of the support.
- > Check the external surface of the cover by using a light source from inside in order to check for holes and/or tears once a year.

### Cleaning and disinfection

- > Between each patient.
- > Wash the cover in a machine at a temperature of up to 90°C.
- > Do not iron the cover.
- > Can be used with a bleach solution of up to 5,000 ppm.

## Information regarding suitable procedures for re-use

- › Each time that it is re-used, the product must be physically and bacteriologically cleaned by performing the procedures given in the technical manual (which can be downloaded at [www.winnicare.com](http://www.winnicare.com)).
- › Replace the cover if the appearance of its surface changes.
- › Do not allow wet surfaces to come into direct contact with the foam.
- › Do not allow scouring, stripping, or “sharps” products to come into direct contact with the protection.

## 5 WARNINGS AND PRECAUTIONS OF USE

### Precautions

If there is a risk of pressure sores, use relief devices or positioning systems or replace the mattress by a therapeutic mattress for the prevention of pressure sores.

You must:

- › Change positions at least once every 4 hours.
- › Maintain skin hygiene and avoid maceration.
- › In case of incontinence, regularly change any protective items.
- › Once a day, observe the condition of the skin, or ensure that it is observed.
- › Ensure that the patient's food intake is sufficient and appropriate.
- › Regularly drink enough fluids.
- › Limit large distances between the body and the support by as much as possible, with the exception of the bed support sheet, clothing, or any adult nappies.
- › Ensure that there are no foreign bodies such as: tubing, crumbs, grease, etc.

### Circumstances under which a user should consult with a health professional

If one of the required measures cannot be followed, you must inform your doctor or nurse as soon as possible. Inform your doctor or nurse as soon as possible of any abnormal events such as fever, pain, reddening or whitening of the weight-bearing areas (head, shoulder, back, hip, shoulder blade, pelvis, hip etc.)

### Warnings

- › Not compatible with external cardiac massage, unless a board is placed between the thorax and the upper surface of the mattress,
- › The 'Medical bed - accessories - non-therapeutic mattress' assembly shall comply with the requirements for the protection against patient entrapment in the moveable parts as per EN 60601-1-52
- › A minimal distance of 22 cm between the top of the bed rail and the top of the mattress with the patient on top of it is advised..
- › This mattress must be installed on a medical bed base which is in good condition and compliant with safety and performance requirements.
- › Enhanced monitoring for patients who are at risk of falls (patient profile: invalid, agitated, not lucid).

### Important

Consult with your prescriber within one month of acquiring this support.

## Lifespan information

The device can be reused as long as its initial performance is maintained throughout its life cycle. Its estimated lifespan is 6 years.

### Accessories

Only the mattress protection accessory supplied by the manufacturer ensures the safety and performance of the complete product.

## 6 STORAGE/ HANDLING/ DISPOSAL

### Specific conditions for storage and/or handling

- › Store away from light, heat, and humidity.
- › The device's rapid identification label, located on the packaging of the device, must always be visible during storage.

### Product disposal

Do not throw the product out into non-dedicated locations in nature. Use the recycling facilities available in your country.

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Mangar International Limited (WINNOCARE group)
Cubro Limited. 149 Taurikura Drive, Tauriko,