
# APLOT & PHARMA PLOTS

APLOT 40 < 120 kg  
88 < 265 lb 90 cm

PHARMA PLOTS 40 < 90 kg  
88 < 198 lb 90, 120, 140 cm

**WINNCARE FRANCE, Miniac**  
ZONE ACTIPOLE  
2 RUE DE SAINT-COULBAN  
35540 MINIAC-MORVAN - FRANCE

**WINNCARE FRANCE, Nîmes**  
200 rue Charles Tellier  
Actiparc de Grézan  
30034 NIMES Cedex 1 - FRANCE

IFU\_WAFM

IFU\_WAFM

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## 1. INSTRUCTIONS

### Purpose of the device

Prevention of pressure ulcers. Handicap compensation.

### Indications

The waffle foam mattresses are bed supports which help to prevent pressure ulcers in patients bedridden for 10 to 15 hours a day, who mobilise themselves with difficulty, with no major neurological disorders, no arterial disease, and in good to average general condition.

- > One or more stage 1 and/or 2 pressure ulcers, not located on the weight bearing areas, or with the possibility of support exclusion (using a localized pressure relief system), patient standing or not during the day.
- > One stage 3 or 4 pressure ulcer not located on the weight bearing areas, with the possibility of support exclusion using a localized pressure relief system,
- > One stage 1 or 2 pressure ulcer in the weight bearing areas, patient standing or not during the day, with a positioning system and intervention by medical assistants 3 times per day for verification of the installation and for turning.

### Contra-indications

- > Periodic repositioning is not possible,
- > Patients at high to very high of pressure ulcers,
- > Patients with stage 3 or 4 pressure ulcers in the weight bearing areas,
- > Patients whose weight is higher than that which the mattress can support.
- > Use on a medical imaging table.

### Patient and user target groups

Adult patients who are hospitalized, institutionalized, or at home, presenting with one or more pressure ulcers and/or who are at risk of developing pressure sores due to any temporary or permanent alterations in their health. These individuals are cared for by health professionals, with the help of assistants if required.


### Indicate the undesirable side-effects

All serious incidents which are related to the device must be notified to the manufacturer and the competent authority of the member state in which the user and/or patient resides.

Inform the competent authority if you believe, or if you have reason to believe, that the device presents a serious risk, or that it has been tampered with.

## 2. COMPOSITION OF THE MEDICAL DEVICE

- 1 Polyurethane foam mattress ( $\geq 34 \text{ kg/m}^3$ ).
- 2 Removable cover with zip (PROMUST), or DERMALON fitted sheet.
- 3 Identification label on the outside of the cover.
- 4 Instructions for use (paper or electronic).

## 3. CLINICAL BENEFITS, PERFORMANCE, MECHANISM OF ACTION

### Performance characteristics of the device

The waffle foam mattresses allow for pressure relief in supporting areas to improve blood circulation and thus reduce the risk of pressure ulcers in the weight bearing areas of the body, Reduces shearing and maceration.

### Expected clinical benefits

Maintenance of tissue oxygenation in the anatomical areas which are in contact with the supporting surface, by reducing the pressure applied to the cutaneous and sub-cutaneous tissue.

### Information for health professionals

Observe the condition of the patient's skin which is in contact with the support structure several times a day. Verify that the mattress is correctly installed.

## 4. USAGE INSTRUCTIONS

### User training

User training must be provided by people who are trained and validated by the relevant economic operators. The choice of support must take into account the risk of pressure ulcers, as per the Braden scale, or any other validated scale, the individual's weight, and its use by the caregivers.

### Installation of the device

When the mattress is supplied with a PROMUST cover, it is delivered zipped. When the mattress is delivered with a DERMALON fitted sheet, the sheet is placed on the mattress, as it is delivered next to it. The surface of the foam made up of studs must be positioned in contact with the patient's body. The protection must be dry. Place the mattress on the bed in the right direction, i.e. with the flap of the PROMUST cover facing downwards or with the sides of the DERMALON fitted sheet/ underpad facing downwards. Lay the patient down.

## MATTRESS DELIVERED COMPRESSED

Foot/head reversible mattress

⌚ 30 min → 2h \*

\* This time may be increased with the storage period (0 to 6 months).

The return to a normal visual appearance of the mattress under normal usage temperature conditions may take up to 24 hours, but this does not present any risk for his use.

## MATTRESS DELIVERED FLAT

Foot/head reversible mattress

## Preventive maintenance

- › Annually carry out a visual inspection of the state of the foam: visible subsidence of the material, and a non-equal and very slow return to the foam's original state are signs of ageing which compromise its pressure ulcers prevention properties.
- › Once a year, check the external surface of the cover, by using a light source from inside in order to check for holes and/or tears.

## Cleaning and disinfection

- › Between each patient.
- › Machine wash the PROMUST cover at up to 90°C. Do not machine wash the DERMALON fitted sheet.
- › Do not iron the cover.
- › Compatible with a bleach solution of up to 5000 ppm (PROMUST), or up to 1000 ppm (DERMALON).

## Information regarding suitable procedures for re-use

- › Each time that it is re-used, the product must be physically and bacteriologically cleaned by performing the procedures given in the technical manual (which can be downloaded at [www.winnicare.com](http://www.winnicare.com)).
- › Replace the cover if the appearance of its surface changes.

## 5. WARNINGS AND PRECAUTIONS

### Required measures

A support by itself is not enough to prevent pressure ulcers. Other preventative measures are also required.

- › Change positions at least once every 4 hours.
- › Maintain skin hygiene and avoid maceration.
- › In case of incontinence, regularly change any protective items.
- › Once a day, observe the condition of the skin, or ensure that it is observed.
- › Ensure that the patient's food intake is sufficient and appropriate.
- › Regularly drink enough fluids.
- › A minimal distance of 22 cm between the top of the bed rail and the top of the mattress with the patient on it is advised.
- › Limit large distances between the body and the support by as much as possible, with the exception of the bed support sheet, clothing, or adult diapers.
- › Ensure that there are no foreign bodies such as: tubing, crumbs, grease, etc.
- › Reinforced surveillance in subjects at risk of falling (patient profile: invalid, agitated, not lucid).

### Circumstances under which a user should consult with a health professional

If one of the required measures cannot be followed, you must inform your doctor or nurse as soon as possible. Inform your doctor or nurse as soon as possible of any abnormal events such as fever, pain, reddening or whitening of the weight-bearing areas (head, shoulder, back, hip, shoulder blade, pelvis, hip etc.)

## Precautions

- › Bone and/or muscle trauma which has not been stabilized and which is in contact with the support.
- › First days following surgery for a pressure ulcer (skin graft or flap).
- › Patient being followed-up at home, without the possibility of receiving help from medical assistants.
- › Also verify the condition of the exposed part of the skin each time that care is provided, and during positional changes.

## Warnings

- › Not compatible with external cardiac massage, unless a board is placed between the thorax and the upper surface of the mattress.
- › Re-evaluate the patient's risk of entrapment in the non-moving parts of the medical bed associated with "therapeutic mattress and accessories and articulated positions of the bed frame" in accordance with the EN 60601-2-52/A1 standard for adults.
- › The use of dry hot water bottles heated in the microwave in direct contact with the surface of the mattress is prohibited.

## Lifespan information

The device can be reused as long as its initial performance is maintained throughout its life cycle. Its estimated lifespan is 6 years.

## Accessories

Only the mattress protection accessory supplied by the manufacturer ensures the safety and performance of the complete product.

## 6. STORAGE/HANDLING/ DISPOSAL

### Specific conditions for storage and/or handling

- › Store away from light, heat, and humidity.

### Disposal of the product

Do not throw the product out into non-dedicated locations in nature. Use the recycling facilities available in your country.

**Pressure (max) in mmHg before ageing, using the maximum validated weight**

	Max. user weight	Max. user pressure		
		Head	Buttocks	Heel
APLIT	120 kg	54	46	92
PHARMA PLOTS	90 kg	65	54	76

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