



## UNIVERSAL AMPUTEE with head support

For users who have a single or double amputation

### USAGE

- For transfers from a sitting or lying position.
- For transfers in mobile and stationary hoists.
- The version with head support is applied for users, who need extra support around the head.

### FEATURES

- The sling is transferring user from a wheelchair/chair or from lying position from a bed or bath.
- The sling provides good support of the body and stabilizes the pelvis, which give the user a feeling of stability and safety.
- Comes with an extra soft cuff which can be fastened around the amputated leg.
- The sling has shared leg straps with special strap solution, which means that the leg straps are not crossed. The leg straps have a large padded surface, which reduces the pressure on the entire femur.
- The sling ensures an upright sitting position that contributes to a correct positioning of the user.
- The sling has a small seat opening.
- Comes with a location pocket for better sling fitting when sitting in a wheelchair, and guide straps on both sides of the sling to make the transfer easier.



### SPECIFICATIONS

SWL: 275 kg

**Material:** Polyester

**Sizes:** XS-S-M-L-XL-XXL

### ITEM NUMBERS:

Universal Amputee, with extra leg strap and head support, polyester:


	Item no.
<b>XS</b>	009-09753
<b>S</b>	009-09754
<b>M</b>	009-09755
<b>L</b>	009-09756
<b>XL</b>	009-09757
<b>XXL</b>	009-09758

## SLING SIZE GUIDE

The sling has a coloured strap on upper back with size and midline mark, which ensures correct placement of the sling.

Colour code is presented on the size chart.

The image below shows a yellow colour code, which means that the sling is in size M (Medium):



	XS	S	M	L	XL	XXL
<b>Y</b> (hip measurement)	37 cm	42 cm	48 cm	52 cm	68 cm	78 cm
<b>Z</b> (measurement from head to lower back)	75 cm	80 cm	91 cm	100 cm	100 cm	100 cm



Please consult the technical manual for further specifications. All figures and measurements are for guidance only and may differ in actual situations. Winncare Nordic reserves the right to changes without notice. Product photos may include accessories. Please consult your dealer for further information.

## FOR MORE INFORMATION VISIT

[winncare.com](http://winncare.com)

E-mail: [info@winncare.dk](mailto:info@winncare.dk)

Phone: +45 70273720

WIN  CARE