

A new comprehensive range with a **PLUS**, suitable for the majority of patients across multiple care settings

17 shapes available dependent on patient need

| Forme | Dénomination | Size in cm (L x D x H or Diam) | |
|-------|---------------------------------|-----------------------------------|-----------------|
| | | CAREWAVE PLUS | POZ'N'FORM PLUS |
| 1 | Universal Cushion | 45 x 30 x 10 | 45 x 30 x 15 |
| 2 | | 55 x 40 x 12 | 55 x 40 x 15 |
| 3 | | 45 x 20 x 8 | 45 x 20 x 10 |
| 4 | Circular Cushion | 50 x 45 x D10 | 50 x 45 x D15 |
| 5 | Semi-circular Cushion | 60 x 35 x D10 | 60 x 35 x D15 |
| 6 | Cylindrical Cushion | 60 x D20 | |
| 7 | | 45 x D12 | |
| 8 | Slanted Cushion | 40 x 20 x 9 | 40 x 20 x 12 |
| 9 | Lateral 30° Positioning Cushion | 180 x 55 x 15 | |
| 10 | | 170 x 55 x 15 | |
| 11 | Half-moon Cushion | 130 x 80 x 20 | |
| 12 | Triangular Cushion | 60 x 30 x 30 | |
| 13 | Small Triangular Cushion | 45 x 25 x 25 | |
| 14 | Semi-fowler Cushion | 70 x 50 x 10 | |
| 15 | Heel off-loading Cushion | 80 x 20 x D15 | |
| 16 | Occipital Rest Cushion | - | 50 x 40 x D10 |
| 17 | Thigh Abduction Cushion | 38 x 30 x 12 | 38 x 30 x 16 |

Choose from 2 types of technology based on clinical need

POZ 'IN' FORM **PLUS**



microfibrés


For **POSITIONING**



Reduced mobility


> **COMFORT**

CAREWAVE **PLUS**



microbeads

For **POSTURE**



Reduced mobility to immobile

> **CARE**

Choose from 3 cover materials bringing together comfort, hygiene and durability to meet the demands of any care settings

Cover **Standard**



REUSABLE COVER

High density polyurethane

> Easy to clean

Cover **Comfort**



REUSABLE COVER

Natural cellulose fibres

> A feeling of comfort and security

Cover **Hygiene**



Polyurethane treated with silver ions ultrasonically welded

> Safe and easy to use for multiple patients



All our materials have received the strictest level, class 1 OEKO-TEX certification, indicating no harmful substances, suitable for the most sensitive skin.

