



ALOVA™ PEDIA CONFORT



Braden scale score



KEY FEATURES

- ▶ Single-piece mattress made up of:
 - Molded visco elastic top layer: density: $\geq 75 \text{ Kg/m}^3$; hardness: 1.8 kPa
 - High Resilient (HR) foam base: density: 28 Kg/m^3 ; hardness: 3.6 kPa
- ▶ Patient weight: 15 to 60 Kg (VMA/PEDIA/S), 20 to 90 Kg (VMA/PEDIA/M)
- ▶ Dimensions: 160 x 80 x 14 cm (62.9 x 31.4 x 5.5 inches)
- ▶ Mattress weight: 12 Kg
- ▶ Cover available: Promust CIC

KEY BENEFITS

- ▶ Visco elastic foam provides excellent pressure re-distribution properties: homogeneous immersion of the body into the foam, thus increasing the bearing surface and reducing pressure applied in areas at risk of pressure ulcers.
- ▶ 2 versions available: for children from 3 to 5 years old (VMA/PEDIA/S) or from 6 to 12 years old (VMA/PEDIA/M).
- ▶ Visco elastic foam provides patient comfort.

INDICATIONS

HELP IN PREVENTING PRESSURE ULCERS FOR PATIENT AT MEDIUM TO HIGH RISK OF PRESSURE ULCER DEVELOPMENT, BEDRIDDEN MORE THAN 15 HOURS A DAY.

HELP IN TREATING EXISTING PRESSURE ULCERS FROM STADE 1 TO 2 COMBINED WITH TECHNICAL AID DEVICES.

PATIENT WEIGHT: 15 TO 90 KG/ 2.3 TO 14.1 ST

- ▶ Single-piece mattress made up of visco elastic foam top layer on High Resilient (HR) polyurethane foam for pediatric medical beds
- ▶ Removable and washable integral cover

DESCRIPTION	ITEM REFERENCES ACCORDING TO PACKAGING CARDBOX	FIRE STANDARDS
ALOVA™ FOR CHILDREN FROM 3 TO 5 YEARS OLD	VMA/PEDIA/S	EN 597-1 and 2 GPEM D1 90 Classe D
ALOVA™ FOR CHILDREN FROM 6 TO 12 YEARS OLD	VMA/PEDIA/M	

3 years warranty

Cleaning instructions: Promust CIC cover



Photos not contractually binding

FP000094/01_2018

